



Scan here for access to customer resources



## Review and apply these helpful, practical tips when using the cobas® Liat® PCR System



### Stable surface reminder

Ensure analyzer is on a surface free of vibration from instruments like a vortex or centrifuge.



### Refrigerate assay tubes

Keep assay tubes chilled until test is ready to be performed



### Sample volume


Use the volume specified in the PI:  
1 mL for Strep (liquid amies)  
3 mL for Flu and Flu & RSV (UTMs)

## Sample collection and handling

Sample collection is an important part of the testing process. Limit contaminants like blood and mucus that can inhibit PCR.<sup>1</sup>

For best results, collect samples within 3-7 days of symptom onset.<sup>2</sup>

### Flu and Flu/RSV testing



Have patients blow nose before nasopharyngeal sampling

Leave swab in nostril for a few seconds to absorb secretions<sup>2</sup>

Swirl swab in media before breaking swab shaft

Viability of pathogens with freeze/thaw diminishes and may lead to false negatives<sup>2</sup>

### Strep A testing



Avoid blood and phlegm

Touch only tonsils and back of throat when collecting sample

Swirl swab in media 3x before breaking swab shaft

Vigorously shake the media tube for 5 seconds to evenly disperse patient specimen

### Pipetting



Keep it sterile

Squeeze pipette bulb fully before dipping in sample & release slowly

Put pipette tip down to first segment of tube (pull up as you release)

Make sure there are no bubbles

### Clean analyzer frequently

Wipe down with soft, lint-free cloth and 70% isopropanol or 5-10% bleach solution.\*



- ✓ Keep product literature readily available
- ✓ Use only transport media from PI

If you have any questions or need assistance with using the cobas® Liat® PCR System, please visit [go.roche.com/liatimplementation](https://go.roche.com/liatimplementation) or call 1 (800) 800-5973

<sup>1</sup> Schrader, C., Schielke, A., Ellerbroek, L. and Johne, R. (2012), PCR inhibitors – occurrence, properties and removal. *J Appl Microbiol*, 113: 1014-1026. Accessed April 24, 2018.  
<sup>2</sup> Influenza Specimen Collection. Center for Disease Control and Prevention. <https://www.cdc.gov/flu/pdf/freeresources/healthcare/flu-specimen-collection-guide.pdf>. Accessed May 7, 2018.)

