

# Deliver GC Workflow Solutions



## Sample Preparation

From sample prep to separation, Fisher Chemical™ GC Resolv™ solvents and solvent blends have set the standard of excellence for consistent, reproducible performance in GC.

Product	Cat. No.
Methanol (GC Resolv), Fisher Chemical	A457-4
Hexanes (GC Resolv), Fisher Chemical, 4L	H307-4
Acetone (GC Resolv), Fisher Chemical, 4L	A928-4
Methylene Chloride (Optima), Fisher Chemical, 4L	D151-4
Toluene (Optima), Fisher Chemical, 4L	T291-4

## Sample Handling

Thermo Scientific™ MS Certified Vials are tested and certified for up to 15 critical physical characteristics. They are pre-cleaned, low particle and low background, leading to precise results and unmatched consistency.

Product	Cat. No.
9mm MS Certified Clear Screw Thread Kits	03-377-300
9mm MS Certified Amber Screw Thread Kits	03-377-301
9mm MS Certified Clear Screw Thread Kits	03-377-306
9mm MS Certified Amber Screw Thread Kits	03-377-307
9mm MS Certified Clear Screw Thread Kits	03-060-050
9mm MS Certified Clear Screw Thread Kits	03-060-051
9mm Certified Clear Wide Opening Screw Thread Kits	03-060-052
9mm Certified Clear Wide Opening Screw Thread Kits	03-060-053

## Separation

Thermo Scientific™ TraceGOLD Fast GC Columns are ID columns that are shorter and smaller than conventional GC columns, meaning analysis times can be reduced, yielding faster analysis with the same separation.

Product	Cat. No.
TraceGold TG-5MS	26-098-1420
TraceGold TG-5MS	26-098-1540
TraceGold TG-Wax MS	26-088-1420
TraceGold TG-1MS	26-099-1420
TraceGold TG-624	26-085-3390
TraceGold TG-5SIIIMS	26-096-1420
TraceGOLD TG-624SIIIMS	26-085-3390

## Analysis

30% reduction in run time over a standard GC Column with no changes in system configuration.

## Questions about Chemicals? Looking for Something Specific?

Our chemical specialists are experts on our entire industry-leading selection of chemicals, and they're ready to answer any specific questions you may have or help you place an order. Just fill out the form on this page and a specialist will contact you.

Visit [fishersci.com/Chromatography](http://fishersci.com/Chromatography) or [fishersci.ca/Chromatography](http://fishersci.ca/Chromatography) for more information.

### In the United States:

For customer service, call 1-800-765-7000  
To fax an order, use 1-800-926-1166  
To order online: [fishersci.com](http://fishersci.com)

### In Canada:

For customer service, call 1-800-234-7437  
To fax an order, use 1-800-463-2996  
To order online: [fishersci.com](http://fishersci.com)

