


Disposable Isolation Gown

- Laminated fabric
- Polypropylene: 28gsm
- Polyethylene: 15gsm
- Knitted cuffs
- Blue
- PPE Category I
- Non-sterile
- Non-surgical
- Back slit
- With neck tie and belt
- Size: Universal



Cat. No	Description	Pack qty
16804471	Disposable Isolation Gown	2000

Dimensions	Specification	Tolerance
Universal		
A. Width, top folded gown	70cm	±2
B. Length, sleeve	62cm	±2
C. Length, gown	118cm	±2
D. Width, arm opening	28cm	±2
E. Length, knitted cuff opening	5,5cm	±1.3
F. Length, neck to bottom	109cm	±2
G. Length, waist tie	175cm	±8.5
H. Length, neck tie	32cm	±2.5
Weight	42gsm	±3

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk
 Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it
 Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl
 Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es
 Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2020 Thermo Fisher Scientific Inc. All rights reserved.
 Trademarks used are owned as indicated at fishersci.com/trademarks.

