



Take it easy— Column packages for ÄKTA™ start

Whether you are purifying antibodies, tagged proteins or untagged (native) proteins we offer a package to meet your needs. Choose from five different start-up application packages, each one consisting of a set of prepacked columns adapted for your purification task.

- User-friendly prepacked columns for shorter start-up time
- Column packages adapted for your purification needs
- High performance media for reliable purification results



Ordering information

Product	Content	Code number*
Tagged package, His <i>Purification of Histidine-tagged proteins</i>	HiTrap™ HP (1 × 1 ml) HiTrap™ TALON® crude (1 × 1 ml) HiTrap Desalting (1 × 5 ml)	29-0588-03
Tagged package, GST <i>Purification of GST-tagged proteins</i>	GSTrap™ 4B (1 × 1 ml) HiTrap Desalting (1 × 5 ml)	29-0588-04
Antibody package, Protein A <i>Purification of human antibodies</i>	HiTrap Protein A HP (1 × 1 ml) HiTrap Desalting (1 × 5 ml)	29-0588-05
Antibody package, Protein G <i>Purification of most species of antibodies, including rat</i>	HiTrap Protein G HP (1 × 1 ml) HiTrap Desalting (1 × 5 ml)	29-0588-06
Untagged package <i>Purification of most proteins, untagged or native</i>	HiTrap Q HP (1 × 1 ml) HiTrap SP HP (1 × 1 ml) HiPrep™ 16/60 Sephacryl™ S-200 HR	29-0588-07

* Code numbers only valid when ordering direct from GE Healthcare. If you are ordering from distributors, please contact your regional sales representative.

Need some guidance in selecting the right package?

Please visit the media selection guide supported on the ÄKTA start homepage:
www.gelifesciences.com/aktastart

For other purification needs the Purify selection guide app will guide you:
www.gelifesciences.com/purify



GE, imagination at work, and GE monogram are trademarks of General Electric Company. ÄKTA, HiTrap, HisTrap, GSTrap, HiPrep, and Sephacryl are trademarks of GE Healthcare companies. TALON is a registered trademark of Clontech Inc. © 2013 General Electric Company - All rights reserved. First published Sep. 2013 GE Healthcare Bio-Sciences AB, Björkgatan 30, 751 84 Uppsala, Sweden