



**Fisherbrand**<sup>®</sup>  
QUALITY. RELIABILITY. VALUE.

## Fisherbrand Carts

*Designed with your lab in mind.*

### Overview

Fisherbrand Laboratory Carts fit in every lab. These versatile carts have endless uses from transporting equipment and supplies to providing auxiliary bench space.

Catalog No.	L/W/H (in)	Description	Wt (lbs)
15252400	34.8 x 19 x 36.4	<b>Stainless Steel Cart.</b> Features two SST pans that are 2" deep on three sides with one open side for easy transfer of equipment from countertop to cart	56
15252401	34.8 x 19 x 36.4	<b>2-Basket Glassware Cart.</b> Used to carry glassware and small instruments in 2 vinyl-coated wire baskets. Includes removable plastic drip pan to catch spills	43
15252402	34.8 x 19 x 36.4	<b>4-Basket Glassware Cart.</b> Used to carry glassware and small instruments in 4 vinyl-coated wire baskets. Includes removable plastic drip pan to catch spills	43
15252403	34.8 x 19 x 36.4	<b>Pan Cart.</b> Used for transport of wet or soiled labware. Features 2" deep powder coated steel, removable watertight pans lined with gray rubber mats to prevent shifting	58
15252404	34.8 x 19 x 36.4	<b>Portable Table.</b> Features two open shelves that provide easy transfer of equipment from countertop to cart	44
15252405	34.8 x 19 x 36.4	<b>Chemical Cart.</b> Features two open shelves with rigid PVC extrusions on two sides to prevent containers from sliding	46

Fisherbrand is a registered trademark of Fisher Scientific.

### All Carts Feature:

- 4" diameter casters, two with toe locks
- Ready-to-roll, full factory assembly
- Ergonomic, extruded rubber handle grips
- Durable Fisherbrand blue, powder-coated 18 gauge welded tubular steel frame
- ID tag
- 5 year warranty
- Supports loads up to 400 lbs.



15252401



15252402



15252400



15252403



15252404



15252405



**Fisher Scientific**

For customer service, call 1-800-766-7000.  
To fax an order, use 1-800-926-1166.  
To order online: [www.fishersci.com](http://www.fishersci.com)

© 2013 Thermo Fisher Scientific Inc. All rights reserved.

BN0103133