

MicroAmp plastic consumables compatibility chart for Applied Biosystems endpoint PCR systems and genetic analyzers

Category/part description	Cat. No.	3 x 32-well	96-well		96-well Fast	384-well		Genetic analyzers		
		ProFlex™	Veriti™, ProFlex, SimpliAmp™	2720	9700	Veriti	ProFlex, Veriti	9700	310	3130, 3130xl, 3500, 3500xL, 3730, 3730xl
96-well 0.2 mL plates										
Optical 96-Well Plate	N8010560, 4316813		•	•	•					•
Optical 96-Well Plate with Barcode	4306737, 4326659		•	•	•					•
96-Well Plate with Barcode & Optical Caps	403012		•	•	•					
Optical 96-Well Plate with Barcode & Optical Adhesive Films	4314320		•	•	•					
EnduraPlate Optical 96-Well Clear Plate with Barcode*	4483354, 4483352		•	•	•					•
96-well 0.1 mL plates										
Fast Optical 96-Well Plate, 0.1 mL	4346907					•			•	•
Fast Optical 96-Well Plate with Barcode, 0.1 mL	4346906, 4366932					•			•	•
EnduraPlate Optical 96-Well Fast Clear Plate with Barcode*	4483485, 4483494					•			•	•
384-well plates										
Optical 384-Well Plate	4343370						•	•		•
Optical 384-Well Plate with Barcode	4309849, 4326270, 4343814						•	•		•
EnduraPlate Optical 384-Well Clear Plate with Barcode*	4483285, 4483273						•	•		•
Strip tubes and caps										
Fast 8-Tube Strip, 0.1 mL	4358293					•				
Optical 8-Tube Strip with attached optical caps, 0.2 mL	A30588	•	•	•	•					
8-Tube Strip with attached domed caps, 0.2 mL	A30589	•	•	•	•					
8-Tube Strip, 0.2 mL	N8010580, N8010838	•	•	•	•				•	
Optical 8-Tube Strip, 0.2 mL	4316567	•	•	•	•					
8-Cap Strip	N8010535, N8011535, N8010835	•	•	•	•	•				
Optical 8-Cap Strip	4323032	•	•	•	•	•				
12-Cap Strip	N8010534, N8011534, N8010834		•	•	•	•				
Single tubes										
Fast Reaction Tube with Cap, 0.1 mL	4358297, 4358293					•				
Reaction Tube with Cap, 0.2 mL	N8010540, N8010612, N8010840, N8011540	•	•	•	•					
Reaction Tube without Cap, 0.2 mL	N8010533, N8011533, N8010833	•	•	•	•					
Optical Tube without Cap, 0.2 mL	N8010933	•	•	•	•					
Seals and covers										
Clear Adhesive Film	4306311		•		•	•	•	•		
Optical Adhesive Film	4360954, 4311971		•		•	•	•	•		
96-Well Full Plate Cover	N8010550			•	•					
Accessories										
Splash Free 96-Well Base	4312063		•	•	•	•				
96-Well Support Base	4379590		•	•	•	•				•
96-Well Base	N8010531		•	•	•					
96-Well Reaction Tube/Tray/Retainer Set, 0.2 mL	403083, 403086			•	•					

*Multiple colors are available.

MicroAmp plastics consumables compatibility chart for Applied Biosystems real-time PCR systems

Category/part description	Cat. No.	48-well	96-well			96-well Fast			384-well
		StepOne™	7000	7300; 7500	QuantStudio™, ViiA7™, 7900HT	StepOne™ Plus	7500	QuantStudio, ViiA™ 7, 7900HT	QuantStudio, ViiA 7, 7900HT
96-well 0.2 mL plates									
Optical 96-Well Plate	N8010560, 4316813		•	•	•				
Optical 96-Well Plate with Barcode	4306737, 4326659		•	•	•				
Optical 96-Well Plate with Barcode & Optical Caps	403012		•	•	•				
Optical 96-Well Plate with Barcode & Optical Adhesive Films	4314320		•	•	•				
EnduraPlate Optical 96-Well Clear Plate with Barcode*	4483354, 4483352		••		•				
96-well 0.1 mL plates									
Fast Optical 96-Well Plate, 0.1 mL	4346907					•	•	•	
Fast Optical 96-Well Plate with Barcode, 0.1 mL	4346906, 4366932					•	•	•	
EnduraPlate Optical 96-Well Fast Clear Plate with Barcode*	4483485, 4483494					•	•	•	
384-well plates									
Optical 384-Well Plate	4343370								•
Optical 384-Well Plate with Barcode	4309849, 4326270, 4343814								•
EnduraPlate Optical 384-Well Clear Plate with Barcode*	4483285, 4483273								•
48-well plates									
Fast Optical 48-Well Plate	4375816	•							
Strip tubes and caps									
Fast 8-Tube Strip, 0.1 mL	4358293	•				•	•	•	
Optical 8-Tube Strip with attached optical caps, 0.2 mL	A30588		•	•	•				
Optical 8-Tube Strip, 0.2 mL	4316567		•	•	•				
Optical 8-Cap Strip	4323032	•	•	•	•	•	•	•	
Single tubes and caps									
Fast Reaction Tube with Cap, 0.1 mL	4358297	•				•			
Optical Tube without Cap, 0.2 mL	N8010933		•	•					
Seals and covers									
Optical Adhesive Film	4360954, 4311971		•	•	•	•	•	•	•
48-Well Optical Adhesive Film	4375323	•							
Reaction trays									
96-Well Tray/Retainer Set	403081		•						
Fast 48-Well Tray	4375282	•							
96-Well Tray for VeriFlex Blocks	4379983					•			
Accessories									
Splash-Free 96-Well Base—10 bases	4312063		•	•	•	•	•	•	
96-Well Support Base—10 bases	4379590		•	•	•	•	•	•	
96-Well Base	N8010531		•	•	•	•	•	•	

*Multiple colors are available. ••• = Requires use of proper adapter, Cat. No. A24820.

Find out more at thermofisher.com/plastics



In the United States:

For customer service, call 1-800-766-7000
 To fax an order, use 1-800-926-1166
 To order online: fishersci.com

In Canada:

For customer service, call 1-800-234-7437
 To fax an order, use 1-800-463-2996
 To order online: fishersci.ca

