



RICCA CHEMICAL COMPANY is the premier choice for High Purity Water. We offer water to meet USP Purified, USP/EP Purified, ASTM, ACS, HPLC and other requirements. Our water is specially purified by conventional water softening, pre-filtration, activated carbon organic adsorption, reverse osmosis, mixed bed double-deionization, ultraviolet light irradiation, and 0.2 micron filtration. Custom options in batch size and packaging are also available.

ACS/ASTM For general laboratory requirements

- Continual production monitoring
- Meets current ACS/ASTM requirements
- ASTM Types I-III

ACS Reagent Grade, ASTM Type I, ASTM Type II

Note: Although this water is not distilled, it has the chemical purity equal to or better than double distilled water.

1L Part Number: 915032

Sizes: 4L, 10L, 20L, 55 Gal

Packaging: Poly Natural, Cubitainer, Ropak, Poly Drum

ACS Reagent Grade, ASTM Type I, ASTM Type II

Note: Packaged in glass containers for applications requiring extremely low levels of organics.

4L Part Number: 91511

Sizes: 4L

Packaging: Glass Amber

ACS Reagent Grade, ASTM Type II

1L Part Number: 915232

Sizes: 500mL, 4L, 10L, 20L, 55 Gal

Packaging: Poly Natural, Cubitainer, Ropak, Poly Drum

ACS Reagent Grade, ASTM III

1L Part Number: 915332

Also available in 4L

Packaging: Glass Amber

HPLC/ACS For general laboratories with LC/HPLC/UPLC testing

- Suitable for Ultraviolet Spectrophotometry and Liquid Chromatography
- Absorbance and LC Suitability results on CofA

HPLC Grade, ACS Reagent Grade

1L Part Number: 915332

Also available in 4L

Packaging: Glass Amber

DISTILLED For users that require distilled water

- Additional distillation step to further purify High Purity Water
- Filled into specially cleaned glass containers

Distilled, Reagent Grade

Note: The high purity water feedstock for this product is manufactured by a process validated to produce water meeting ACS and ASTM specifications.

1L Part Number: 918032

Sizes: 4L

Packaging: Glass Amber

USP/EP PURIFIED For pharmaceutical testing labs

- Water system validated under cGMP
- Tested to the latest USP/EP
- Low microbial/Non-sterile

USP/EP Purified

4L Part Number: 91901

Sizes: 10L, 20L

Packaging: Poly Natural, Glass Amber, Cubitainer

USP Purified, Non-Sterile

Note: There is no microbial specification for USP Purified Water

Part Number: 919100055E

Packaging: Poly Drum



In the United States:

For customer service, call 1-800-766-7000.

To fax an order, use 1-800-926-1166.

To order online: www.fishersci.com

In Canada:

For customer service, call 1-800-234-7437.

To fax an order, use 1-800-463-2996.

To order online: www.fishersci.ca

RIGHT • READY • RICCA